

Trendelenburg Ergonomic chair

AJ18 Chair

Part number: A091802



 Add: 1590 S Milliken Ave., unit A, Ontario, CA 91761

 www.adsdental.com www.adsequip.com

 E-mail: sales@adsequip.com

 Tel: 626-620-0456

 Technical Support: 800-488-9708 (West Coast: 626-633-7636 East Coast: 626-366-3848)

 Facebook: ADS Dental

 Youtube: ADS Dental System

 Twitter: ADS Dental System

 Instagram: ADS Dental

Standard Features

Trendelenburg ergonomic chair,
AJ18 hydraulic chair with 3 programmable positions,
Double articulating headrest,
Swing up/down iron armrest,
Cast iron chair frame,
Iron backrest stand,


Cast aluminum base,
Wide backrest, upgraded fiber leather,
60 degrees seat rotation, chair safety switch,
One touchpad, Foot control



Part number

A091802

Upholstery Color Selections



Specifications

Weight	340lbs
Max patient load	550lbs
Electrical requirements	AC110V±22V; 50HZ,60HZ;
Environmental requirements	Ambient temperature 41°F~104°F; Humidity≤80%; Atmospherical pressure 8.7psi~15.4psi
Motor run time	Intermittent operation 2minOn / 18min Off

Dimensions

Unit: inch

