

OSO Pure's SC-Series Competitive Set

Criteria	SC-100 OSO Pure	SC-400 OSO Pure	VistaPure Crosstex	Sterisil AC+ Solmetex	Notes
Retail price	\$3,239	\$5,997	\$3,554	\$3,612	
Pressure pump (integrated)	Y	Y	\$500*	\$500*	Overcomes low tap water pressure, maximum TDS removal & refresh rate. *Price for pump, tank pressure sensor & LOPSI solenoid, installation.
Refresh rate (gal/hr)	4	16	Up to 4	Up to 2	"Up to" indicates best possible results should local water conditions conform to manufacturer's specifications of pressure & quality.
Tank Capacity (gal)	3.2 ¹	6.0 ¹	4.25	2.8	¹ SC Series high refresh rates support instrument washer with smaller tank, easier to place in cabinetry. Larger tanks available from all.
Quick-change filters	Y	Y	N	Y	Easy filter change by staff.
Filter change interval	1 year ¹	1 year ¹	Varies by filter ²	Varies by filter ²	¹ Except for practices with very high TDS tap water, change all SC-100 filters once annually (SCHEDULED). ² Other systems exhaust filters randomly (DISRUPTIONS) depending on water quality & usage.
Complete PM kit (2026)	\$440	\$654	\$703 ¹	\$877 ¹	¹ 2024 pricing
Bottle fill faucet	Included	Included	Included	Included	
Autoclave fill wand	Included	Included	Included	Included	
Instrument washers	1	4+	1 ¹	0 ²	¹ S3010 fittings kit required ² Requires tank upgrade + fittings
Phone support	Factory techs	Factory techs	Call center script	Call center script	

Call OSO Pure @ 800-676-7873