



DESIGNED TO FIT. BUILT TO PERFORM.

WWW.BOWWOOD.NET



At Bow Wood Dental Systems, we believe your operator cabinetry should maximize efficiency and enhance your clinic.

The operator is the focal point of every practice, so ensure all your patients feel welcome and cared for.

We offer an 88" tall option for our Rear Support Cabinets for increased privacy in an open concept area and for extra storage.



Our Rear Support unit is offered in two depths, 19" or 23". This practical option allows you to accommodate wheelchairs and comply with handicap regulations in smaller than standard size operator rooms.

Rear Support Treatment Centre



PULLOUT SURFACE
Strong. Stable. Versatile. Those are just three words to describe our pullout work surface. Push it in. It disappears. Pull it out and your assistant has an extra work area that locks. Ask for our dual pullout option.



DUO SWING TABLE
Did you know that 1 in 10 dentists are left-handed? Take advantage of our Duo Swing table to include ambidextrous capability at every chair, or simply utilize the full range of chairside positioning it offers.



TIP-OUT BINS
Customize your Rear Support Treatment center with quick access bins that open to a 45° angle. Locate items in seconds and check inventory at a glance.



Dual L-R Sliding Counters



Single Pullout



Linear Sliding Table

Centre Island Console



X-RAY ENCLOSURE

Conceal the x-ray machine and have ops private at the same time. The back end of our centre island can include the necessary cabinetry to do both. A popular choice.



CENTRE ISLAND COLUMN

Our laser cut steel ladder frame construction gives exceptional rigidity so X-rays don't drift, and bolt patterns are compatible with lights from any manufacturer.



STORAGE OPTION


The island rear section can be specified as a storage area configured to your needs. A capacious solution to an ongoing problem.



Centre Island Round



Upper & Base Storage



X-Ray Pass Thru



Maximize your office space by using Centre Island Consoles in your office design. You can share equipment and a sink between rooms. And save space by removing the wall.

Bow Wood's Slimline Centre Console is just 22" wide, saving you even more space, and improving your office workflow.



Side Treatment Centre



ROUND END

We round our end corner to give you more room to get in and out of your operatory and improve wheelchair accessibility . Plus if we do say so ourselves, it looks really smart.



GLOVE/CUP/TOWEL DISPENSER

Our GCT upper storage cabinets come with time-saving dispensers so you never have to open a door to get what you need and risk contamination.



PERFECT DRAWER


Customize the drawer sizes to accommodate specific equipment such as an ultrasonic scaler unit or electric handpiece.




Weave purpose and beauty together to create a system that works for your practice.

Give your dental office a trendy new look by using vessel sinks. With a lowered countertop so the sinks remains at standard height, it's user friendly, and it looks great.


Or, if there's no need for an extra sink, our side cabinets can be customized to remove the sink to save on plumbing costs.




Glove/Cup/Towel Dispenser



Storage Upper - 1 door



Storage Upper - 2 door



Storage Upper - 3 door

The health and safety of your patients is paramount in your practice, so effective and efficient processing of instruments in the sterilization centre is essential.

The 6 main steps of the Instrument Processing Protocol (cleaning, inspecting, packaging, sterilization, sterile storage and quality assurance) happen within five modules of the sterilization centre.

Through intuitive design and process management you can create a system of fully custom dental cabinets that complement your workflow perfectly.



Sterilization Centre



MAXIMUM TRAY STORAGE

We designed wire racks so they give the maximum amount of storage space possible. They're chrome, long-lasting... and they're adjustable.




BLUE OR RED LIGHTS

Our sterilization centres can also be outfitted with blue and red lights for increased patient awareness of the aseptic process. So, give us the green light to put blue and red lights in for you.



SOLID SURFACE PULLOUTS

Solidly waterproof. Our tower pullouts make sure any moisture that drips down during sterilization stays exactly where it's supposed to. Because containment is where it's at.



SC 12

SC L-Shape

SC Galley Style

Mobile Cart



MCA-18



MCU-18



MCO



Alabama



North-Carolina

Create flexible storage spaces, a movable workspace and increase the efficiency of your operatory rooms by using mobile carts. Bow Wood's mobile carts come in multiple sizes and configuration options to create a solution that is perfect for your space. Our ortho carts are specially designed to suit orthodontic and pedodontics practices and support those procedures and storage requirements. They can be used as a fully functional treatment station with the integrated delivery system and storage options. Can be either mobile or fixed depending on your practice needs.





Serving all of North America

www.bowwood.net

design@bowwood.net

1-888-269-9663



February 2024—Printed in Canada