



LEO



LEO 8'/9'/10' Ceiling mount LED Dental Surgical Light

Part number		Retail Price
A0601651	8' Ceiling mount	\$4,499
A0601652	9' Ceiling mount	\$4,499
A0601653	10' Ceiling mount	\$4,499

LEO Mobile LED Dental Surgical Light

Part number	Retail Price
A0605651	\$4,499

Specifications


- Illuminance: 30,000-120000lux
- Color rendering index: $\geq 92Ra$
- Color temperature: $4000 \pm 500k$
- Spot diameter: $3.94in \leq d_{10} \leq 9.84in$
- Length of light: $\geq 70in$
- Ceiling static load: 330 lbs
- Floor height: $\geq 9.52ft$
- Light height adjustment range: 38.2in-86.6in

- Illuminance: 30,000-120000lux
- Color rendering index: $\geq 92Ra$
- Color temperature: $4000 \pm 500k$
- Spot diameter: $3.94in \leq d_{10} \leq 9.84in$
- Length of light: $\geq 70in$
- Floor height: $\geq 9.52ft$
- Light height adjustment range: 29.5in-77.0in


Dimensions

Unit: inch

LEO Ceiling mount LED Dental
Surgical Light



LEO Mobile LED Dental
Surgical Light





LIBRA



Libra 8'/9'/10' Ceiling mount LED Dental Light

Part number		Retail Price
A0601641	8' Ceiling mount	\$2,699
A0601642	9' Ceiling mount	\$2,699
A0601643	10' Ceiling mount	\$2,699

Libra Mobile LED Dental Light

Part number	Retail Price
A0605641	\$2,699

Specifications


- Illuminance: 17,500-70000lux
- Color rendering index: $\geq 90Ra$
- Color temperature: $4000 \pm 500k$
- Spot diameter: $3.94in \leq d_{10} \leq 9.84in$
- Length of light: $\geq 70in$
- Ceiling static load: 140lbs
- Floor height: $\geq 9.52ft$
- Light height adjustment range: 55in-106in

- Illuminance: 17,500-70000lux
- Color rendering index: $\geq 90Ra$
- Color temperature: $4000 \pm 500k$
- Spot diameter: $3.94in \leq d_{10} \leq 9.84in$
- Length of light: $\geq 70in$
- Floor height: $\geq 9.52ft$
- Light height adjustment range: 36in-86in

Dimensions

Unit: inch

Libra Ceiling mount
LED Dental Light



The size of "L" depend on the height of the room

Libra Mobile LED
Dental Light

